

## **BOAR Autocorrelator**

BOAR stands for Bimirror based Optical Autocorrelation with Retrieval. This new technique of ultrashort pulses characterization relies on interferometric single shot autocorrelation and two photon absorption. The time delay is encoded into a spatial interferogram which is used to evaluate the pulse duration, the  $2\omega$  spectrum and the chirp. There is no non-linear crystal and no phase matching issues, the spectral working range is therefore very broad (1200 - 2400 nm). The BOAR is actually combining all the advantages : simple, extremely robust, accurate and reliable measurements, spatially resolved, suitable for rather chirped pulses and the retrieval is done in real time directly by Fourier transformation. Two models are available for two different temporal windows.



- Very easy to use
- Temporal and spectral measurement
- Real time chirp measurement
- No phase matching issues
- Broad spectral range¹
- Nonsensitive to polarization
- Suitable for any rep rate
- Single shot up to 150 kHz<sup>2</sup>
- Sub-10 fs in the NIR
- Achromatic and non-dispersive

Models	BOAR FS	BOAR PS	1. with a single optics set.
Pulse duration range	Sub-10 to 750 fs	Sub-10 to 1500 fs	The spectral range is directly accessible and there is no need for any manipulations.
Temporal window Δτ	3 ps	5.5 ps	
Spectral range (nm)	1200 - 2400 nm¹		2. If the laser rep rate is higher than the single shot capacity of the BOAR, the
Shot to shot measurement capacity	150 kHz with synchronisation 75 kHz without²	30 kHz with synchronisation 15 kHz without²	measurements are average over several shots. For example when measuring a 1MHz laser with the BOAR
Input pulse repetition rate	From Hz to GHz		FS, it is not possible to measure less than 12 pulses together. Devices
Input pulse energy and average power (for 100 fs pulses)	Single sho 1 MHz : 5	Standard models Single shot: 100 nJ 1 MHz: 5 nJ / 5 mW 100 MHz: 200 pJ / 20 mW with higher single measurement capacit be made upon reques	
	With low energy option Single shot : 5 nJ 1 MHz : 100 pJ / 100 μW 100 MHz : 10 pJ / 1 mW 1 GHz : 1 pJ / 1 mW		
Input polarization	any		0.5 (1.5 (1.0 2010
Detection	CMOS 12 bit - 6 Mpx - 72 dB - USB 3.1	CMOS 12 bit - 18 Mpx - 72 dB - USB 3.1	© Femto Easy 09-2018.  Product specifications are subject to change
Dimensions (mm)	125×150×80		without prior notice.